SpherHA

SpherHA Dense Granules

SpherHA Porous Chips

SpherHA Injectable Paste

SpherHA Mouldable Crunch

	Non usare se la confezione è danneggiata Do not use if package is damaged No lo use si el paquete está dañado Nie używaj, jeśli opakowanie jest uszkodzone
	Non resterilizzare Do not resterilize No volver a esterilizar Nie sterylizuj ponownie

IT

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Eumed S.r.l.
Via Tommasco, 68 - 35131 Padova (ITALIA)

EN

SpherHA - line of synthetic bone substitutes based on nano-structured bio-mimetic hydroxyapatite

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SpherHA - Sterile, single-use device

SpherHA Dense Granules

SpherHA Porous Chips

SpherHA Injectable Paste

SpherHA Mouldable Crunch

INTENDED USE

SpherHA is a line of synthetic bone substitutes based on nano-structured bio-mimetic hydroxyapatite for applications in dental and oral and maxillo-facial surgery, orthopaedics, traumatology and neurosurgery.

DESCRIPTION

SpherHA hydroxyapatite is a calcium phosphate compound equivalent to the mineral matrix of human bone, both in chemical composition and in structure and size of its crystals, with a Ca/P ratio of 1.67 similar to that of natural bone apatite. The porous and interconnected structure of the hydroxyapatite composing SpherHA devices has an optimal osteoconductive activity, promoting cell penetration and colonization, circulation of nutrients and rapid vascularization. In addition, their nano-structure is characterized by a large surface area in relation to volume, such as to promote the adhesion of the cells that invade it, also allowing its complete degradation by osteoclasts and a complete remodeling into new vital bone tissue.

AVAILABLE FORMATS

The devices in the SpherHA line are available in the form of dense granules and porous chips, in both cases with a granulometry of 0.5-1 mm and 1-2 mm. They are also available in injectable paste and mouldable crunch.

PRODUCT CONSTITUENTS

SpherHA Dense Granules

Dense granules of nano-structured bio-mimetic hydroxyapatite with granulometry of 0.5-1 mm and 1-2 mm.
SpherHA Porous Chips
Porous chips of nano-structured bio-mimetic hydroxyapatite with granulometry of 0.5-1 mm and 1-2 mm.
SpherHA Injectable Paste
Injectable paste with high density of nano-crystals and micro-powder of nano-structured bio-mimetic hydroxyapatite in a buffered aqueous solution at physiological pH.
SpherHA Mouldable Crunch
High-density mouldable crunch of nano-crystals, micro-powder and micro-granules of nano-structured bio-mimetic hydroxyapatite in a buffered aqueous solution at physiological pH.

INDICATIONS

Bone substitutes of the SpherHA line are indicated for bone loss reconstruction procedures and more specifically in the filling of bone defects of traumatic, degenerative, pathological or surgical origin, in association, where required, with osteosynthesis means and/or prosthetic implants. Dense granules and porous chips can also be mixed or associated with the following biological fluids or autologous tissues: peripheral blood, platelet rich plasma (PRP), bone marrow blood, bone marrow aspirated cell concentrates, bone tissue or adipose tissue. SpherHA may also be used in combination with hemostats, fibrin glue, growth factors and antibiotics.

CONTRAINDICATIONS AND SIDE EFFECTS

SpherHA is a biocompatible device. Latex and Phthalates free. SpherHA should not be used in patients with proven hypersensitivity to calcium phosphate salts. It is also contraindicated in the presence of inflammation or infection at the implant site or in adjacent tissues.

INSTRUCTIONS FOR USE

The following instructions in no way replace the usual surgical procedures for the treatment of bone defects. The SpherHA line of bone substitutes, in all their formats, are ready to use. Once the box is opened, remove the envelope containing the product. Open it in aseptic way, either by handing the contents to sterile personnel or by dropping it in a sterile field. Dense granules and porous chips can be mixed before use, as indicated in the instructions.

Prepare the receiving site taking care to remove any fibrous and/or necrotic tissue and bleeding bone surfaces. Then apply the product to the implant site in an adequate quantity to properly fill the defect. Gently compact the porous granules with beaters or other dedicated instruments. Remove any excess material to avoid dislocation of the graft. When using SpherHA in injectable paste or mouldable crunch format, do not wash in the surgical area after implantation to avoid dispersion and/or runoff. In oral surgery, it is advisable to cover the grafts with guided bone regeneration membranes (Collygen Membrane / Collygen RS) to avoid fibrotic tissue infiltration at the implant site. The SpherHA line of bone substitutes does not have adequate mechanical strength for implantation in sites subject to loading and in these applications must therefore always be associated with suitable metal fixation or osteosynthesis. Degradation and remodeling, in absence of external factors that contribute to accelerating or delaying the process, usually occurs between 6-8 months for SpherHA Porous Chips, Injectable Paste and Mouldable Crunch, and 12-14 months for SpherHA Dense Granules.

WARNINGS AND PRECAUTIONS

Read the instructions carefully before the use. SpherHA is a Class III, single-use, single-patient medical device and cannot be reused or re-sterilized. Any unused product residue must not be stored and used in other operations or for other patients and must be disposed of as contaminated hospital waste in accordance with hospital regulations. In the event that the disposable device is mistakenly reused, as well as affecting its performance, possible infections could occur to the patient due to the lost sterility and/or possible contamination of the device itself. Failure to comply with this indication will result in a different use than that intended by the manufacturer, who can no longer guarantee the conditions of hygiene, functionality and safety once the packaging has been opened.

Bone substitutes of the SpherHA line do not have mechanical characteristics suitable to withstand direct loads and must, in these cases, always be associated with fixation systems or metal osteosynthesis. SpherHA has not been tested on pregnant patients. Do not use the device if the package is open or damaged. Do not use SpherHA after the expiration date. The expiration date refers to the product in its undamaged and regularly stored packaging. Keep out of the reach of children.

PRODUCT PACKAGING

SpherHA is supplied STERILE.

Dense granules

Bottle made of polyethylene terephthalate (PET), polyethylene (PE), in laminated polyester-aluminium (ALU/OPA-PE) envelope, single or double. Illustrative leaflet.

Porous chips

Bottle in polyethylene terephthalate (PET), polyethylene (PE), in ALU/OPA-PE envelope, single or double. Illustrative leaflet.
Syringe made of polypropylene (PP), polyethylene (PE) with "Filter Cap" (screw filter cap, removable), in ALU/OPA-PE envelope, single. Illustrative leaflet.

Injectable paste

Syringe made of polyamide (PA), polybutylene (PBT), polypropylene (PP), polyethylene (PE), co-polyester (CPE), polycarbonate (PC), with luer-lock connection in ALU/OPA-PE envelope, single or double. Illustrative leaflet.

The package includes a tip for dental formats and a 10 cm long-cannula for formats dedicated to other surgical specialties.

Mouldable crunch

Syringe made of polyamide (PA), polybutylene (PBT), polypropylene (PP), polyethylene (PE), open-mouth screw and luer-lock connection in ALU/OPA-PE envelope, single or double. Illustrative leaflet.

For all formats: the outer label of the product is placed on the outer box (final packaging), while on the outside of each envelope contained in the box there are 4 detachable adhesive patient labels made for traceability, to be applied on the patient record or other documentation.
Envelope breakage and disposal of the package
Do not use the product if the package is damaged. The package materials do not require special disposal conditions.

STERILIZATION

SpherHA bone substitutes are sterilized using 25 kGy gamma rays.

STORAGE

The device must be stored in a cool and dry place at a temperature between 4 and 40°C avoiding direct exposure to sunlight. If stored properly, the integrity of the packaging and the sterility of the device are guaranteed for 3 years from the date of production (see expiration date on the outer label).

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Tiss’You S.r.l.