



Product photo

|                       |  |
|-----------------------|--|
| Manufacturing company | <ul style="list-style-type: none"> <li>▪ Clover Orthopedics s.r.l.</li> </ul>  |
| Distributor company   | <ul style="list-style-type: none"> <li>▪ Clover Orthopedics s.r.l.</li> </ul>  |
| Trading name          | <ul style="list-style-type: none"> <li>▪ Rally Percutaneous Drills</li> </ul>  |
| Product type          | Rally is a device intended to be used for the removal of bone tissue in the field of foot surgery  |
| Materials             | <ul style="list-style-type: none"> <li>▪ Stainless steel 17- 4-Ph (AISI 630)</li> </ul>  |
| Indications           | <p>Treatment, through minimally invasive percutaneous access, of the metatarsophalangeal joint using cutting tools (burs).</p> <p>The application of the devices in question must be carried out by medical-surgeon staff with excellent knowledge of orthopedic surgery for foot joint pathologies.<br/>The choice of the use of surgical instruments and possible and practicable treatments is exclusive to the medical surgeon.<br/>Any surgical decisions other than those recommended by the manufacturer are at the discretion and responsibility of the surgeon.</p>   |
| Contraindications     | <p>Its use is not recommended if present in the patient:</p> <ul style="list-style-type: none"> <li>▪ severe local and/or systemic infections;</li> <li>▪ allergies to metals such as stainless steel or nickel;</li> <li>▪ pathologies of the circulatory system;</li> <li>▪ excess weight;</li> <li>▪ diabetes.</li> </ul> <p>If the Clover Orthopedics Rally Burs are considered the best solution for the patient, but the patient has one or more of the aforementioned contraindications, it is essential that the patient is informed of the negative consequences that could hinder the success of the surgical procedure. For more information, see the instructions for use of the device in question.<br/>Le frese Rally sono disponibili nelle seguenti tipologie:</p> |

## Product features

- Shannon isham short and standard
- Wedge
- Cylindrical in helix
- Straight cylindrical
- Conical
- Brophy

available in different diameters and cutting lengths.

The drills, selected by the surgeon, must be made available and removed from their sterile packaging shortly before their use.

It is always necessary to check that they are firmly connected to the rotating motor spindle, making sure before use that they are blocked and not deformed. The manufacturer recommends manual use.

The cutters must never exceed 8000 revolutions per minute. The drills, during their use, should be cooled by means of physiological water by gravity.

After the surgery, make sure that the patient follows an adequate post-operative course, paying particular attention to the hygiene of the surgical accesses performed.

### *Notes for Disposal*

Used devices must be disposed of in accordance with the procedures for potentially contaminated material.

The package must not be dispersed in the environment.

## Packaging and sterilization

Each package consists of a double Tyvek-PE bag and an external box.  
RALLY burs are for single use only.

All the devices, in the different variants, are supplied STERILE. They undergo a sterilization process using ethylene oxide. If the package is damaged, it is recommended not to use its contents.

It is not intended that the devices can be re-sterilized.

## Brand

- Clover Orthopedics

## Class

- IIa

## Certifications

- CE Mark (0426) | Certificate 517-00-00-DM



| CND   | RDM       | CODICE         | DESCRIZIONE   |
|---|-----------|----------------|---|
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025686/R | RLY-A00SS00122 | Fresa Rally SHANNON ISHAM<br>d 1.2 L 70 mm          |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025698/R | RLY-A00SS10133 | Fresa Rally WEDGE<br>d 3.1 L 70 mm                  |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025699/R | RLY-A10SS10134 | Fresa Rally WEDGE<br>d 4.1 L 70 mm                  |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025700/R | RLY-A00SS08082 | Fresa Rally SHANNON ISHAM CORTA<br>d 1.2 L 70 mm    |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025701/R | RLY-A00SS15122 | Fresa Rally CILINDRICA IN ELICA"<br>d 2 L 70 mm     |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025702/R | RLY-A10SS10133 | Fresa Rally CONICAL BURR<br>d 3.1 L 70 mm           |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025703/R | RLY-A00SS13155 | Fresa Rally BROPHY<br>d 5 L 70 mm                   |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025704/R | RLY-A00SS00202 | Fresa Rally CYLIN.L STRAIGHT BURR<br>d 2 L 70 mm    |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2097815/R | RLY-A00SS00152 | Fresa Rally CYLIN.L STRAIGHT BURR<br>d 2x15 L 70 mm |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025705/R | RLY-B00SS00122 | Fresa Rally SHANNON ISHAM<br>d 1.2 L 65 mm          |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025706/R | RLY-B00SS10133 | Fresa Rally WEDGE<br>d 3.1 L 65 mm                  |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025707/R | RLY-B10SS10134 | Fresa Rally WEDGE<br>d 4.1 L 65 mm                  |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025708/R | RLY-B00SS08082 | Fresa Rally SHANNON ISHAM CORTA<br>d 1.2 L 65 mm    |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025709/R | RLY-B00SS15122 | Fresa Rally CILINDRICA IN ELICA<br>d 2 L 65 mm      |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025710/R | RLY-B10SS10133 | Fresa Rally CONICAL BURR<br>d 3.1 L 65 mm           |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025711/R | RLY-B00SS13155 | Fresa Rally BROPHY<br>d 5 L 65 mm                   |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025712/R | RLY-B00SS00202 | Fresa Rally CYLIN.L STRAIGHT BURR<br>d 2 L 65 mm    |

| CND   | RDM       | CODICE         | DESCRIZIONE   |
|---|-----------|----------------|---|
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2097873/R | RLY-B00SS00152 | Fresa Rally CYLIN.L STRAIGHT BURR<br>d 2x15 L 65 mm |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025713/R | RLY-C00SS00122 | Fresa Rally SHANNON ISHAM<br>d 1.2 L 60 mm          |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025714/R | RLY-C00SS10133 | Fresa Rally WEDGE<br>d 3.1 L 60 mm                  |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025715/R | RLY-C10SS10134 | Fresa Rally WEDGE<br>d 4.1 L 60 mm                  |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025716/R | RLY-C00SS08082 | Fresa Rally SHANNON ISHAM CORTA<br>d 1.2 L 60 mm    |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025717/R | RLY-C00SS15122 | Fresa Rally CILINDRICA IN ELICA<br>d 2 L 60 mm      |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025718/R | RLY-C10SS10133 | Fresa Rally CONICAL BURR<br>d 3.1 L 60 mm           |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025719/R | RLY-C00SS13155 | Fresa Rally BROPHY<br>d 5 L 60 mm                   |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025720/R | RLY-C00SS00202 | Fresa Rally CYLIN.L STRAIGHT BURR<br>d 2 L 60 mm    |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2097880/R | RLY-C00SS00152 | Fresa Rally CYLIN.L STRAIGHT BURR<br>d 2x15 L 60 mm |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025721/R | RLY-D00SS00122 | Fresa Rally SHANNON ISHAM<br>d 1.2 L 55 mm          |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025722/R | RLY-D00SS10133 | Fresa Rally WEDGE<br>d 3.1 L 55 mm                  |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025723/R | RLY-D10SS10134 | Fresa Rally WEDGE<br>d 4.1 L 55 mm                  |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025724/R | RLY-D00SS08082 | Fresa Rally SHANNON ISHAM CORTA<br>d 1.2 L 55 mm    |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025725/R | RLY-D00SS15122 | Fresa Rally CILINDRICA IN ELICA<br>d 2 L 55 mm      |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025726/R | RLY-D10SS10133 | Fresa Rally CONICAL BURR<br>d 3.1 L 55 mm           |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025727/R | RLY-D00SS13155 | Fresa Rally BROPHY<br>d 5 L 55 mm                   |

| CND   | RDM       | CODICE         | DESCRIZIONE   |
|---|-----------|----------------|---|
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2025728/R | RLY-D00SS00202 | Fresa Rally CYLIN.L STRAIGHT BURR<br>d 2 L 55 mm    |
| P091301 - FRESE MONOUSO<br>PER PROTESICA ORTOPEDICA | 2097885/R | RLY-D00SS00152 | Fresa Rally CYLIN.L STRAIGHT BURR<br>d 2x15 L 55 mm |