



Product photo

- Manufacturing company ▪ Clover Orthopedics s.r.l.
- Distributor company ▪ Clover Orthopedics s.r.l.
- Trading name ▪ Misano

Product type The Misano thoracolumbosacral stabilisation system from Clover Orthopedics consists of single- and polyaxial pedicle screws, clamping nuts, straight and pre-curved bars, laminar and pedicle hooks, and cross-links and is used to promote solid arthrodesis of the part of the spine being treated.

- Materials**
- Implantable titanium alloy (Ti6Al4V) conforming to ASTM F-136 and ISO 5832-3
 - Implantable Titanium conforming to ASTM F-67 and ISO 5832-2
 - Chrome Cobalt alloy conforming to ASTM F-1537 (contained exclusively in MISANO Cr-Co Bar)

Indications Properly used, the Misano thoracolumbar-sacral stabilisation system from Clover Orthopedics is indicated to promote the development of a solid thoracic, lumbar and sacral arthrodesis. It is recommended in cases of spinal deformities, degenerative disc disease, traumatic vertebral fractures, vertebral tumours, spinal stenosis, spondylolisthesis, pseudoarthrosis and previous unsuccessful attempts at vertebral arthrodesis. Any surgical decisions other than those recommended by the manufacturer are at the discretion and responsibility of the surgeon. Do not use 4.5 mm diameter screws in the lumbar and lumbosacral spine and do not couple 4.5 mm diameter screws with Cr-Co rods.

Contraindications Contraindications to the implantation of Clover Orthopedics' Misano posterior spinal stabilisation system are similar to those of similar products on the market, and include, but are not limited to:

Absolute contraindications

- Active infections;
- Possible allergies to the metals indicated in the normative references in the "Materials" section;
- Allergy to Nickel, Chromium, Cobalt (when using MISANO Cr-Co Bar);

- Uncooperative patients unable to follow instructions.

For further details on the chemical compositions of the Misano system materials please refer to the regulations listed in this data sheet

Relative contraindications

- Metastases;
- Severe muscular, neurological, or vascular disease;
- Fever or leukocytosis;
- Pregnancy, except for treatment of unstable vertebral fractures;
- Signs of inflammation at the implant site;
- Inadequate soft tissue coverage at the site of surgery;
- Severe osteoporosis.

If the implantation of the Misano system by Clover Orthopedics is considered to be the best solution for the patient, and if the patient has one or more of the contraindications listed above, it is essential to inform the patient of the possible negative consequences for the success of the operation.

Product features

- The System consists of polyaxial and uniaxial screws with a progressive thread, i.e. with a wider pitch in the distal part of the screw, in order to speed up its entry into the vertebral body, and shorter in the proximal part to improve pedicle anchorage.
- The polyaxial screw of the Misano system is designed to allow both open (OPEN) and minimally invasive (MIS) surgical approaches, being able to connect to dedicated cannulas that allow percutaneous insertion.
- Misano system screws are available both solid and cannulated and fenestrated to allow the injection of bone cement and improve the primary stability of the system in the case of osteoporotic bone.
- The Misano system instrumentation includes dedicated cannulas that connect to the head of the screw and allow it to be inserted percutaneously. The cannulas also allow, in open surgical approaches, the use of dedicated persuaders that allow a fine and precise recovery of vertebral misalignments such as spondylolisthesis.
- The MSN spinal stabilisation system from Clover Orthopedics is intended to be Nickel-Free with the exception of the cobalt chrome bars (MISANO Cr-Co Bar).

Packaging and sterilization

The system is washed and packaged inside a sealed double medical PE+Tyvek pouch suitable for ETO sterilisation.

Brand

Clover Orthopedics

Class

IIb

Certifications

CE Mark (0426)

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015865	MSN-A00T500000	MISANO - Cap di serraggio per viti e uncini
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015238	MSN-B00T545025	MISANO - Vite monoassiale - l.25mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015418/R	MSN-B00T545030	MISANO - Vite monoassiale - l.30mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015419/R	MSN-B00T545035	MISANO - Vite monoassiale - l.35mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015420/R	MSN-B00T545040	MISANO - Vite monoassiale - l.40mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015421/R	MSN-B00T545045	MISANO - Vite monoassiale - l.45mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015422/R	MSN-B00T545050	MISANO - Vite monoassiale - l.50mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015423/R	MSN-B00T545055	MISANO - Vite monoassiale - l.55mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015424/R	MSN-B00T555030	MISANO - Vite monoassiale - l.30mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015425/R	MSN-B00T555035	MISANO - Vite monoassiale - l.35mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015426/R	MSN-B00T555040	MISANO - Vite monoassiale - l.40mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015427/R	MSN-B00T555045	MISANO - Vite monoassiale - l.45mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015428/R	MSN-B00T555050	MISANO - Vite monoassiale - l.50mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015429/R	MSN-B00T555055	MISANO - Vite monoassiale - l.55mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015430/R	MSN-B00T565030	MISANO - Vite monoassiale - l.30mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015431/R	MSN-B00T565035	MISANO - Vite monoassiale - l.35mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015432/R	MSN-B00T565040	MISANO - Vite monoassiale - l.40mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015433/R	MSN-B00T565045	MISANO - Vite monoassiale - l.45mm - d 6,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015434/R	MSN-B00T565050	MISANO - Vite monoassiale - l.50mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015435/R	MSN-B00T565055	MISANO - Vite monoassiale - l.55mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015436/R	MSN-B00T575030	MISANO - Vite monoassiale - l.30mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015437/R	MSN-B00T575035	MISANO - Vite monoassiale - l.35mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015438/R	MSN-B00T575040	MISANO - Vite monoassiale - l.40mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015439/R	MSN-B00T575045	MISANO - Vite monoassiale - l.45mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015440/R	MSN-B00T575050	MISANO - Vite monoassiale - l.50mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015441/R	MSN-B00T575055	MISANO - Vite monoassiale - l.55mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015442/R	MSN-B00T585030	MISANO - Vite monoassiale - l.30mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015443/R	MSN-B00T585035	MISANO - Vite monoassiale - l.35mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015444/R	MSN-B00T585040	MISANO - Vite monoassiale - l.40mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015445/R	MSN-B00T585045	MISANO - Vite monoassiale - l.45mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015446/R	MSN-B00T585050	MISANO - Vite monoassiale - l.50mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015447/R	MSN-B00T585055	MISANO - Vite monoassiale - l.55mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015448/R	MSN-B10T545025	MISANO - Vite Monoassiale Cannulata Fenestrata - l.25mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015449/R	MSN-B10T545030	MISANO - Vite Monoassiale Cannulata Fenestrata - l.30mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015450/R	MSN-B10T555030	MISANO - Vite Monoassiale Cannulata Fenestrata - l.30mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015451/R	MSN-B10T565030	MISANO - Vite Monoassiale Cannulata Fenestrata - l.30mm - d 6,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015452/R	MSN-B10T575030	MISANO - Vite Monoassiale Cannulata Fenestrata - l.30mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015453/R	MSN-B10T585030	MISANO - Vite Monoassiale Cannulata Fenestrata - l.30mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015454/R	MSN-B10T545035	MISANO - Vite Monoassiale Cannulata Fenestrata - l.35mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015455/R	MSN-B10T555035	MISANO - Vite Monoassiale Cannulata Fenestrata - l.35mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015456/R	MSN-B10T565035	MISANO - Vite Monoassiale Cannulata Fenestrata - l.35mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015457/R	MSN-B10T575035	MISANO - Vite Monoassiale Cannulata Fenestrata - l.35mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015458/R	MSN-B10T585035	MISANO - Vite Monoassiale Cannulata Fenestrata - l.35mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015459/R	MSN-B10T545040	MISANO - Vite Monoassiale Cannulata Fenestrata - l.40mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015460/R	MSN-B10T555040	MISANO - Vite Monoassiale Cannulata Fenestrata - l.40mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015461/R	MSN-B10T565040	MISANO - Vite Monoassiale Cannulata Fenestrata - l.40mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015462/R	MSN-B10T575040	MISANO - Vite Monoassiale Cannulata Fenestrata - l.40mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015463/R	MSN-B10T585040	MISANO - Vite Monoassiale Cannulata Fenestrata - l.40mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015464/R	MSN-B10T545045	MISANO - Vite Monoassiale Cannulata Fenestrata - l.45mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015465/R	MSN-B10T555045	MISANO - Vite Monoassiale Cannulata Fenestrata - l.45mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015466/R	MSN-B10T565045	MISANO - Vite Monoassiale Cannulata Fenestrata - l.45mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015467/R	MSN-B10T575045	MISANO - Vite Monoassiale Cannulata Fenestrata - l.45mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015468/R	MSN-B10T585045	MISANO - Vite Monoassiale Cannulata Fenestrata - l.45mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015469/R	MSN-B10T545050	MISANO - Vite Monoassiale Cannulata Fenestrata - l.50mm - d 4,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015470/R	MSN-B10T555050	MISANO - Vite Monoassiale Cannulata Fenestrata - l.50mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015471/R	MSN-B10T565050	MISANO - Vite Monoassiale Cannulata Fenestrata - l.50mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015472/R	MSN-B10T575050	MISANO - Vite Monoassiale Cannulata Fenestrata - l.50mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015473/R	MSN-B10T585050	MISANO - Vite Monoassiale Cannulata Fenestrata - l.50mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015474/R	MSN-B10T545055	MISANO - Vite Monoassiale Cannulata Fenestrata - l.55mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015475/R	MSN-B10T555055	MISANO - Vite Monoassiale Cannulata Fenestrata - l.55mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015476/R	MSN-B10T565055	MISANO - Vite Monoassiale Cannulata Fenestrata - l.55mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015477/R	MSN-B10T575055	MISANO - Vite Monoassiale Cannulata Fenestrata - l.55mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015478/R	MSN-B10T585055	MISANO - Vite Monoassiale Cannulata Fenestrata - l.55mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015483/R	MSN-C00T545020	MISANO - Vite Multiassiale - l.20mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015484/R	MSN-C00T545025	MISANO - Vite Multiassiale - l.25mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015485/R	MSN-C00T565030	MISANO - Vite Multiassiale - l.30mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015486/R	MSN-C00T575030	MISANO - Vite Multiassiale - l.30mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015487/R	MSN-C00T585030	MISANO - Vite Multiassiale - l.30mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015488/R	MSN-C00T545030	MISANO - Vite Multiassiale - l.30mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015489/R	MSN-C00T555030	MISANO - Vite Multiassiale - l.30mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015490/R	MSN-C00T565035	MISANO - Vite Multiassiale - l.35mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015491/R	MSN-C00T575035	MISANO - Vite Multiassiale - l.35mm - d 7,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015492/R	MSN-C00T585035	MISANO - Vite Multiassiale - l.35mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015493/R	MSN-C00T595035	MISANO - Vite Multiassiale - l.35mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015494/R	MSN-C00T545035	MISANO - Vite Multiassiale - l.35mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015495/R	MSN-C00T555035	MISANO - Vite Multiassiale - l.35mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015496/R	MSN-C00T565040	MISANO - Vite Multiassiale - l.40mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015497/R	MSN-C00T575040	MISANO - Vite Multiassiale - l.40mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015498/R	MSN-C00T585040	MISANO - Vite Multiassiale - l.40mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015499/R	MSN-C00T595040	MISANO - Vite Multiassiale - l.40mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015500/R	MSN-C00T545040	MISANO - Vite Multiassiale - l.40mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015501/R	MSN-C00T555040	MISANO - Vite Multiassiale - l.40mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015502/R	MSN-C00T565045	MISANO - Vite Multiassiale - l.45mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015503/R	MSN-C00T575045	MISANO - Vite Multiassiale - l.45mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015504/R	MSN-C00T585045	MISANO - Vite Multiassiale - l.45mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015505/R	MSN-C00T595045	MISANO - Vite Multiassiale - l.45mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015506/R	MSN-C00T545045	MISANO - Vite Multiassiale - l.45mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015507/R	MSN-C00T555045	MISANO - Vite Multiassiale - l.45mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015508/R	MSN-C00T565050	MISANO - Vite Multiassiale - l.50mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015509/R	MSN-C00T575050	MISANO - Vite Multiassiale - l.50mm - d 7,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015510/R	MSN-C00T585050	MISANO - Vite Multiassiale - l.50mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015511/R	MSN-C00T595050	MISANO - Vite Multiassiale - l.50mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015512/R	MSN-C00T545050	MISANO - Vite Multiassiale - l.50mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015513/R	MSN-C00T555050	MISANO - Vite Multiassiale - l.50mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015514/R	MSN-C00T565055	MISANO - Vite Multiassiale - l.55mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015515/R	MSN-C00T575055	MISANO - Vite Multiassiale - l.55mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015516/R	MSN-C00T585055	MISANO - Vite Multiassiale - l.55mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015517/R	MSN-C00T595055	MISANO - Vite Multiassiale - l.55mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015518/R	MSN-C00T545055	MISANO - Vite Multiassiale - l.55mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015519/R	MSN-C00T555055	MISANO - Vite Multiassiale - l.55mm - d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015520/R	MSN-C00T565060	MISANO - Vite Multiassiale - l.60mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015521/R	MSN-C00T575060	MISANO - Vite Multiassiale - l.60mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015522/R	MSN-C00T585060	MISANO - Vite Multiassiale - l.60mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015523/R	MSN-C00T595060	MISANO - Vite Multiassiale - l.60mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015524/R	MSN-C00T565065	MISANO - Vite Multiassiale - l.65mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015525/R	MSN-C00T575065	MISANO - Vite Multiassiale - l.65mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015526/R	MSN-C00T585065	MISANO - Vite Multiassiale - l.65mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015527/R	MSN-C00T595065	MISANO - Vite Multiassiale - l.65mm - d 9,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015528/R	MSN-C00T565070	MISANO - Vite Multiassiale - l.70mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015529/R	MSN-C00T575070	MISANO - Vite Multiassiale - l.70mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015530/R	MSN-C00T585070	MISANO - Vite Multiassiale - l.70mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015531/R	MSN-C00T595070	MISANO - Vite Multiassiale - l.70mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015532/R	MSN-C00T565075	MISANO - Vite Multiassiale - l.75mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015533/R	MSN-C00T575075	MISANO - Vite Multiassiale - l.75mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015534/R	MSN-C00T585075	MISANO - Vite Multiassiale - l.75mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015535/R	MSN-C00T595075	MISANO - Vite Multiassiale - l.75mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015536/R	MSN-C00T565080	MISANO - Vite Multiassiale - l.80mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015537/R	MSN-C00T575080	MISANO - Vite Multiassiale - l.80mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015538/R	MSN-C00T585080	MISANO - Vite Multiassiale - l.80mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015539/R	MSN-C00T595080	MISANO - Vite Multiassiale - l.80mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015540/R	MSN-C00T565085	MISANO - Vite Multiassiale - l.85mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015541/R	MSN-C00T575085	MISANO - Vite Multiassiale - l.85mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015542/R	MSN-C00T585085	MISANO - Vite Multiassiale - l.85mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015543/R	MSN-C00T595085	MISANO - Vite Multiassiale - l.85mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015544/R	MSN-C00T565090	MISANO - Vite Multiassiale - l.90mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015545/R	MSN-C00T575090	MISANO - Vite Multiassiale - l.90mm - d 7,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015546/R	MSN-C00T585090	MISANO - Vite Multiassiale - l.90mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015547/R	MSN-C00T595090	MISANO - Vite Multiassiale - l.90mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015548/R	MSN-C10T545020	MISANO - Vite Multiassiale Cannulata Fenestrata - l.20mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015549/R	MSN-C10T545025	MISANO - Vite Multiassiale Cannulata Fenestrata - l.25mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015550/R	MSN-C10T545030	MISANO - Vite Multiassiale Cannulata Fenestrata - l.30mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015551/R	MSN-C10T555030	MISANO - Vite Multiassiale Cannulata Fenestrata - l.30mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015552/R	MSN-C10T565030	MISANO - Vite Multiassiale Cannulata Fenestrata - l.30mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015553/R	MSN-C10T575030	MISANO - Vite Multiassiale Cannulata Fenestrata - l.30mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015554/R	MSN-C10T585030	MISANO - Vite Multiassiale Cannulata Fenestrata - l.30mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015555/R	MSN-C10T545035	MISANO - Vite Multiassiale Cannulata Fenestrata - l.35mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015556/R	MSN-C10T555035	MISANO - Vite Multiassiale Cannulata Fenestrata - l.35mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015557/R	MSN-C10T565035	MISANO - Vite Multiassiale Cannulata Fenestrata - l.35mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015558/R	MSN-C10T575035	MISANO - Vite Multiassiale Cannulata Fenestrata - l.35mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015559/R	MSN-C10T585035	MISANO - Vite Multiassiale Cannulata Fenestrata - l.35mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015560/R	MSN-C10T595035	MISANO - Vite Multiassiale Cannulata Fenestrata - l.35mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015561/R	MSN-C10T545040	MISANO - Vite Multiassiale Cannulata Fenestrata - l.40mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015562/R	MSN-C10T555040	MISANO - Vite Multiassiale Cannulata Fenestrata - l.40mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015563/R	MSN-C10T565040	MISANO - Vite Multiassiale Cannulata Fenestrata - l.40mm- d 6,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015564/R	MSN-C10T575040	MISANO - Vite Multiassiale Cannulata Fenestrata - l.40mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015565/R	MSN-C10T585040	MISANO - Vite Multiassiale Cannulata Fenestrata - l.40mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015566/R	MSN-C10T595040	MISANO - Vite Multiassiale Cannulata Fenestrata - l.40mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015567/R	MSN-C10T545045	MISANO - Vite Multiassiale Cannulata Fenestrata - l.45mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015568/R	MSN-C10T555045	MISANO - Vite Multiassiale Cannulata Fenestrata - l.45mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015569/R	MSN-C10T565045	MISANO - Vite Multiassiale Cannulata Fenestrata - l.45mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015570/R	MSN-C10T575045	MISANO - Vite Multiassiale Cannulata Fenestrata - l.45mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015571/R	MSN-C10T585045	MISANO - Vite Multiassiale Cannulata Fenestrata - l.45mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015572/R	MSN-C10T595045	MISANO - Vite Multiassiale Cannulata Fenestrata - l.45mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015573/R	MSN-C10T545050	MISANO - Vite Multiassiale Cannulata Fenestrata - l.50mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015574/R	MSN-C10T555050	MISANO - Vite Multiassiale Cannulata Fenestrata - l.50mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015575/R	MSN-C10T565050	MISANO - Vite Multiassiale Cannulata Fenestrata - l.50mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015576/R	MSN-C10T575050	MISANO - Vite Multiassiale Cannulata Fenestrata - l.50mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015577/R	MSN-C10T585050	MISANO - Vite Multiassiale Cannulata Fenestrata - l.50mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015578/R	MSN-C10T595050	MISANO - Vite Multiassiale Cannulata Fenestrata - l.50mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015579/R	MSN-C10T545055	MISANO - Vite Multiassiale Cannulata Fenestrata - l.55mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015580/R	MSN-C10T555055	MISANO - Vite Multiassiale Cannulata Fenestrata - l.55mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015581/R	MSN-C10T565055	MISANO - Vite Multiassiale Cannulata Fenestrata - l.55mm- d 6,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015583/R	MSN-C10T575055	MISANO - Vite Multiassiale Cannulata Fenestrata - l.55mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015584/R	MSN-C10T585055	MISANO - Vite Multiassiale Cannulata Fenestrata - l.55mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015585/R	MSN-C10T595055	MISANO - Vite Multiassiale Cannulata Fenestrata - l.55mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015586/R	MSN-C10T565060	MISANO - Vite Multiassiale Cannulata Fenestrata - l.60mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015587/R	MSN-C10T575060	MISANO - Vite Multiassiale Cannulata Fenestrata - l.60mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015588/R	MSN-C10T585060	MISANO - Vite Multiassiale Cannulata Fenestrata - l.60mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015589/R	MSN-C10T595060	MISANO - Vite Multiassiale Cannulata Fenestrata - l.60mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015590/R	MSN-C10T565065	MISANO - Vite Multiassiale Cannulata Fenestrata - l.65mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015591/R	MSN-C10T575065	MISANO - Vite Multiassiale Cannulata Fenestrata - l.65mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015592/R	MSN-C10T585065	MISANO - Vite Multiassiale Cannulata Fenestrata - l.65mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015593/R	MSN-C10T595065	MISANO - Vite Multiassiale Cannulata Fenestrata - l.65mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015594/R	MSN-C10T565070	MISANO - Vite Multiassiale Cannulata Fenestrata - l.70mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015595/R	MSN-C10T575070	MISANO - Vite Multiassiale Cannulata Fenestrata - l.70mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015596/R	MSN-C10T585070	MISANO - Vite Multiassiale Cannulata Fenestrata - l.70mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015597/R	MSN-C10T595070	MISANO - Vite Multiassiale Cannulata Fenestrata - l.70mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015598/R	MSN-C10T565075	MISANO - Vite Multiassiale Cannulata Fenestrata - l.75mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015599/R	MSN-C10T575075	MISANO - Vite Multiassiale Cannulata Fenestrata - l.75mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015600/R	MSN-C10T585075	MISANO - Vite Multiassiale Cannulata Fenestrata - l.75mm- d 8,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015601/R	MSN-C10T595075	MISANO - Vite Multiassiale Cannulata Fenestrata - l.75mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015602/R	MSN-C10T565080	MISANO - Vite Multiassiale Cannulata Fenestrata - l.80mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015603/R	MSN-C10T575080	MISANO - Vite Multiassiale Cannulata Fenestrata - l.80mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015604/R	MSN-C10T585080	MISANO - Vite Multiassiale Cannulata Fenestrata - l.80mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015605/R	MSN-C10T595080	MISANO - Vite Multiassiale Cannulata Fenestrata - l.80mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015606/R	MSN-C10T565085	MISANO - Vite Multiassiale Cannulata Fenestrata - l.85mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015607/R	MSN-C10T575085	MISANO - Vite Multiassiale Cannulata Fenestrata - l.85mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015608/R	MSN-C10T585085	MISANO - Vite Multiassiale Cannulata Fenestrata - l.85mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015609/R	MSN-C10T595085	MISANO - Vite Multiassiale Cannulata Fenestrata - l.85mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015610/R	MSN-C10T565090	MISANO - Vite Multiassiale Cannulata Fenestrata - l.90mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015611/R	MSN-C10T575090	MISANO - Vite Multiassiale Cannulata Fenestrata - l.90mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015612/R	MSN-C10T585090	MISANO - Vite Multiassiale Cannulata Fenestrata - l.90mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015613/R	MSN-C10T595090	MISANO - Vite Multiassiale Cannulata Fenestrata - l.90mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051420/R	MSN-L10T545020	MISANO - Vite Multiassiale Percutanea Fenestrata - l.20mm-d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051421/R	MSN-L10T545025	MISANO - Vite Multiassiale Percutanea Fenestrata - l.25mm-d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051422/R	MSN-L10T565030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051423/R	MSN-L10T575030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051424/R	MSN-L10T585030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm - d 8,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051425/R	MSN-L10T545030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm - d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051426/R	MSN-L10T555030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm-d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051427/R	MSN-L10T565035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051428/R	MSN-L10T575035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051429/R	MSN-L10T585035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051430/R	MSN-L10T595035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051431/R	MSN-L10T545035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm-d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051432/R	MSN-L10T555035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm-d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051433/R	MSN-L10T565040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051434/R	MSN-L10T575040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051435/R	MSN-L10T585040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051436/R	MSN-L10T595040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051437/R	MSN-L10T545040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm-d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051438/R	MSN-L10T555040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm-d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051439/R	MSN-L10T565045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051440/R	MSN-L10T575045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051441/R	MSN-L10T585045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051442/R	MSN-L10T595045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm - d 9,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051443/R	MSN-L10T545045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm-d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051444/R	MSN-L10T555045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm-d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051445/R	MSN-L10T565050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051446/R	MSN-L10T575050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051447/R	MSN-L10T585050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051448/R	MSN-L10T595050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051449/R	MSN-L10T545050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm-d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051450/R	MSN-L10T555050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm-d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051451/R	MSN-L10T565055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051452/R	MSN-L10T575055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051453/R	MSN-L10T585055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051454/R	MSN-L10T595055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051455/R	MSN-L10T545055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm-d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051456/R	MSN-L10T555055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm-d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051457/R	MSN-L10T565060	MISANO - Vite Multiassiale Percutanea Fenestrata - l.60mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051458/R	MSN-L10T575060	MISANO - Vite Multiassiale Percutanea Fenestrata - l.60mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051459/R	MSN-L10T585060	MISANO - Vite Multiassiale Percutanea Fenestrata - l.60mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051460/R	MSN-L10T595060	MISANO - Vite Multiassiale Percutanea Fenestrata - l.60mm - d 9,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051461/R	MSN-L10T565065	MISANO - Vite Multiassiale Percutanea Fenestrata - l.65mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051462/R	MSN-L10T575065	MISANO - Vite Multiassiale Percutanea Fenestrata - l.65mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051463/R	MSN-L10T585065	MISANO - Vite Multiassiale Percutanea Fenestrata - l.65mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051464/R	MSN-L10T595065	MISANO - Vite Multiassiale Percutanea Fenestrata - l.65mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051465/R	MSN-L10T565070	MISANO - Vite Multiassiale Percutanea Fenestrata - l.70mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051466/R	MSN-L10T575070	MISANO - Vite Multiassiale Percutanea Fenestrata - l.70mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051467/R	MSN-L10T585070	MISANO - Vite Multiassiale Percutanea Fenestrata - l.70mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051468/R	MSN-L10T595070	MISANO - Vite Multiassiale Percutanea Fenestrata - l.70mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051469/R	MSN-L10T565075	MISANO - Vite Multiassiale Percutanea Fenestrata - l.75mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051470/R	MSN-L10T575075	MISANO - Vite Multiassiale Percutanea Fenestrata - l.75mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051471/R	MSN-L10T585075	MISANO - Vite Multiassiale Percutanea Fenestrata - l.75mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051472/R	MSN-L10T595075	MISANO - Vite Multiassiale Percutanea Fenestrata - l.75mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051473/R	MSN-L10T565080	MISANO - Vite Multiassiale Percutanea Fenestrata - l.80mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051474/R	MSN-L10T575080	MISANO - Vite Multiassiale Percutanea Fenestrata - l.80mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051475/R	MSN-L10T585080	MISANO - Vite Multiassiale Percutanea Fenestrata - l.80mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051476/R	MSN-L10T595080	MISANO - Vite Multiassiale Percutanea Fenestrata - l.80mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051477/R	MSN-L10T565085	MISANO - Vite Multiassiale Percutanea Fenestrata - l.85mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051478/R	MSN-L10T575085	MISANO - Vite Multiassiale Percutanea Fenestrata - l.85mm - d 7,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051479/R	MSN-L10T585085	MISANO - Vite Multiassiale Percutanea Fenestrata - l.85mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051480/R	MSN-L10T595085	MISANO - Vite Multiassiale Percutanea Fenestrata - l.85mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051481/R	MSN-L10T565090	MISANO - Vite Multiassiale Percutanea Fenestrata - l.90mm - d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051482/R	MSN-L10T575090	MISANO - Vite Multiassiale Percutanea Fenestrata - l.90mm - d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051483/R	MSN-L10T585090	MISANO - Vite Multiassiale Percutanea Fenestrata - l.90mm - d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051484/R	MSN-L10T595090	MISANO - Vite Multiassiale Percutanea Fenestrata - l.90mm - d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051485/R	MSN-L10T545020	MISANO - Vite Multiassiale Percutanea Fenestrata - l.20mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051486/R	MSN-L10T545025	MISANO - Vite Multiassiale Percutanea Fenestrata - l.25mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051487/R	MSN-L10T545030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051488/R	MSN-L10T555030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051489/R	MSN-L10T565030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051490/R	MSN-L10T575030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051491/R	MSN-L10T585030	MISANO - Vite Multiassiale Percutanea Fenestrata - l.30mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051492/R	MSN-L10T545035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051493/R	MSN-L10T555035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051494/R	MSN-L10T565035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051495/R	MSN-L10T575035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051496/R	MSN-L10T585035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm- d 8,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051497/R	MSN-L10T595035	MISANO - Vite Multiassiale Percutanea Fenestrata - l.35mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051498/R	MSN-L10T545040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051499/R	MSN-L10T555040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051500/R	MSN-L10T565040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051501/R	MSN-L10T575040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051502/R	MSN-L10T585040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051503/R	MSN-L10T595040	MISANO - Vite Multiassiale Percutanea Fenestrata - l.40mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051504/R	MSN-L10T545045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051505/R	MSN-L10T555045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051506/R	MSN-L10T565045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051507/R	MSN-L10T575045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051508/R	MSN-L10T585045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051509/R	MSN-L10T595045	MISANO - Vite Multiassiale Percutanea Fenestrata - l.45mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051510/R	MSN-L10T545050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051511/R	MSN-L10T555050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051512/R	MSN-L10T565050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051513/R	MSN-L10T575050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051514/R	MSN-L10T585050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm- d 8,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051515/R	MSN-L10T595050	MISANO - Vite Multiassiale Percutanea Fenestrata - l.50mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051516/R	MSN-L10T545055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm- d 4,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051517/R	MSN-L10T555055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm- d 5,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051518/R	MSN-L10T565055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051519/R	MSN-L10T575055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051520/R	MSN-L10T585055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051521/R	MSN-L10T595055	MISANO - Vite Multiassiale Percutanea Fenestrata - l.55mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051522/R	MSN-L10T565060	MISANO - Vite Multiassiale Percutanea Fenestrata - l.60mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051523/R	MSN-L10T575060	MISANO - Vite Multiassiale Percutanea Fenestrata - l.60mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051524/R	MSN-L10T585060	MISANO - Vite Multiassiale Percutanea Fenestrata - l.60mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051525/R	MSN-L10T595060	MISANO - Vite Multiassiale Percutanea Fenestrata - l.60mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051526/R	MSN-L10T565065	MISANO - Vite Multiassiale Percutanea Fenestrata - l.65mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051527/R	MSN-L10T575065	MISANO - Vite Multiassiale Percutanea Fenestrata - l.65mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051528/R	MSN-L10T585065	MISANO - Vite Multiassiale Percutanea Fenestrata - l.65mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051529/R	MSN-L10T595065	MISANO - Vite Multiassiale Percutanea Fenestrata - l.65mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051530/R	MSN-L10T565070	MISANO - Vite Multiassiale Percutanea Fenestrata - l.70mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051531/R	MSN-L10T575070	MISANO - Vite Multiassiale Percutanea Fenestrata - l.70mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051532/R	MSN-L10T585070	MISANO - Vite Multiassiale Percutanea Fenestrata - l.70mm- d 8,5mm con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051533/R	MSN-L10T595070	MISANO - Vite Multiassiale Percutanea Fenestrata - l.70mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051534/R	MSN-L10T565075	MISANO - Vite Multiassiale Percutanea Fenestrata - l.75mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051535/R	MSN-L10T575075	MISANO - Vite Multiassiale Percutanea Fenestrata - l.75mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051536/R	MSN-L10T585075	MISANO - Vite Multiassiale Percutanea Fenestrata - l.75mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051537/R	MSN-L10T595075	MISANO - Vite Multiassiale Percutanea Fenestrata - l.75mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051538/R	MSN-L10T565080	MISANO - Vite Multiassiale Percutanea Fenestrata - l.80mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051539/R	MSN-L10T575080	MISANO - Vite Multiassiale Percutanea Fenestrata - l.80mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051540/R	MSN-L10T585080	MISANO - Vite Multiassiale Percutanea Fenestrata - l.80mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051541/R	MSN-L10T595080	MISANO - Vite Multiassiale Percutanea Fenestrata - l.80mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051542/R	MSN-L10T565085	MISANO - Vite Multiassiale Percutanea Fenestrata - l.85mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051543/R	MSN-L10T575085	MISANO - Vite Multiassiale Percutanea Fenestrata - l.85mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051544/R	MSN-L10T585085	MISANO - Vite Multiassiale Percutanea Fenestrata - l.85mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051545/R	MSN-L10T595085	MISANO - Vite Multiassiale Percutanea Fenestrata - l.85mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051546/R	MSN-L10T565090	MISANO - Vite Multiassiale Percutanea Fenestrata - l.90mm- d 6,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051547/R	MSN-L10T575090	MISANO - Vite Multiassiale Percutanea Fenestrata - l.90mm- d 7,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051548/R	MSN-L10T585090	MISANO - Vite Multiassiale Percutanea Fenestrata - l.90mm- d 8,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2051549/R	MSN-L10T595090	MISANO - Vite Multiassiale Percutanea Fenestrata - l.90mm- d 9,5mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015614/R	MSN-D00CC00100	MISANO - Barra Cr-Co - d 5.5 mm L.100

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015615/R	MSN-D00CC00105	MISANO - Barra Cr-Co - d 5.5 mm L.105
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015616/R	MSN-D00CC00110	MISANO - Barra Cr-Co - d 5.5 mm L.110
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015617/R	MSN-D00CC00115	MISANO - Barra Cr-Co - d 5.5 mm L.115
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015618/R	MSN-D00CC00120	MISANO - Barra Cr-Co - d 5.5 mm L.120
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015619/R	MSN-D00CC00125	MISANO - Barra Cr-Co - d 5.5 mm L.125
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015620/R	MSN-D00CC00130	MISANO - Barra Cr-Co - d 5.5 mm L.130
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015621/R	MSN-D00CC00135	MISANO - Barra Cr-Co - d 5.5 mm L.135
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015622/R	MSN-D00CC00140	MISANO - Barra Cr-Co - d 5.5 mm L.140
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015623/R	MSN-D00CC00145	MISANO - Barra Cr-Co - d 5.5 mm L.145
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015624/R	MSN-D00CC00150	MISANO - Barra Cr-Co - d 5.5 mm L.150
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015625/R	MSN-D00CC00155	MISANO - Barra Cr-Co - d 5.5 mm L.155
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015626/R	MSN-D00CC00160	MISANO - Barra Cr-Co - d 5.5 mm L.160
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015627/R	MSN-D00CC00165	MISANO - Barra Cr-Co - d 5.5 mm L.165
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015628/R	MSN-D00CC00170	MISANO - Barra Cr-Co - d 5.5 mm L.170
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015629/R	MSN-D00CC00175	MISANO - Barra Cr-Co - d 5.5 mm L.175
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015630/R	MSN-D00CC00180	MISANO - Barra Cr-Co - d 5.5 mm L.180
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015631/R	MSN-D00CC00185	MISANO - Barra Cr-Co - d 5.5 mm L.185
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015632/R	MSN-D00CC00190	MISANO - Barra Cr-Co - d 5.5 mm L.1900

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015633/R	MSN-D00CC00195	MISANO - Barra Cr-Co - d 5.5 mm L.195
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015634/R	MSN-D00CC00200	MISANO - Barra Cr-Co - d 5.5 mm L.200
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015635/R	MSN-D00CC00210	MISANO - Barra Cr-Co - d 5.5 mm L.210
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015636/R	MSN-D00CC00220	MISANO - Barra Cr-Co - d 5.5 mm L.220
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015637/R	MSN-D00CC00230	MISANO - Barra Cr-Co - d 5.5 mm L.230
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015638/R	MSN-D00CC00240	MISANO - Barra Cr-Co - d 5.5 mm L.240
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015639/R	MSN-D00CC00025	MISANO - Barra Cr-Co - d 5.5 mm L.25
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015640/R	MSN-D00CC00250	MISANO - Barra Cr-Co - d 5.5 mm L.250
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015641/R	MSN-D00CC00260	MISANO - Barra Cr-Co - d 5.5 mm L.260
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015642/R	MSN-D00CC00270	MISANO - Barra Cr-Co - d 5.5 mm L.270
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015643/R	MSN-D00CC00280	MISANO - Barra Cr-Co - d 5.5 mm L.280
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015644/R	MSN-D00CC00290	MISANO - Barra Cr-Co - d 5.5 mm L.290
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015645/R	MSN-D00CC00030	MISANO - Barra Cr-Co - d 5.5 mm L.30
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015646/R	MSN-D00CC00300	MISANO - Barra Cr-Co - d 5.5 mm L.300
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015647/R	MSN-D00CC00310	MISANO - Barra Cr-Co - d 5.5 mm L.310
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015648/R	MSN-D00CC00320	MISANO - Barra Cr-Co - d 5.5 mm L.320
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015649/R	MSN-D00CC00330	MISANO - Barra Cr-Co - d 5.5 mm L.330
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015650/R	MSN-D00CC00340	MISANO - Barra Cr-Co - d 5.5 mm L.340

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015651/R	MSN-D00CC00035	MISANO - Barra Cr-Co - d 5.5 mm L.35
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015652/R	MSN-D00CC000350	MISANO - Barra Cr-Co - d 5.5 mm L.350
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015653/R	MSN-D00CC000360	MISANO - Barra Cr-Co - d 5.5 mm L.360
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015654/R	MSN-D00CC000370	MISANO - Barra Cr-Co - d 5.5 mm L.370
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015655/R	MSN-D00CC000380	MISANO - Barra Cr-Co - d 5.5 mm L.380
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015656/R	MSN-D00CC000390	MISANO - Barra Cr-Co - d 5.5 mm L.390
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015657/R	MSN-D00CC000040	MISANO - Barra Cr-Co - d 5.5 mm L.40
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015658/R	MSN-D00CC000400	MISANO - Barra Cr-Co - d 5.5 mm L.400
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015659/R	MSN-D00CC000410	MISANO - Barra Cr-Co - d 5.5 mm L.410
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015660/R	MSN-D00CC000420	MISANO - Barra Cr-Co - d 5.5 mm L.420
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015661/R	MSN-D00CC000430	MISANO - Barra Cr-Co - d 5.5 mm L.430
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015662/R	MSN-D00CC000440	MISANO - Barra Cr-Co - d 5.5 mm L.440
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015663/R	MSN-D00CC000045	MISANO - Barra Cr-Co - d 5.5 mm L.45
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015664/R	MSN-D00CC000450	MISANO - Barra Cr-Co - d 5.5 mm L.450
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015665/R	MSN-D00CC000460	MISANO - Barra Cr-Co - d 5.5 mm L.460
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015666/R	MSN-D00CC000470	MISANO - Barra Cr-Co - d 5.5 mm L.470
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015667/R	MSN-D00CC000480	MISANO - Barra Cr-Co - d 5.5 mm L.480
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015668/R	MSN-D00CC000490	MISANO - Barra Cr-Co - d 5.5 mm L.490

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015669/R	MSN-D00CC00050	MISANO - Barra Cr-Co - d 5.5 mm L.50
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015670/R	MSN-D00CC00050	MISANO - Barra Cr-Co - d 5.5 mm L.500
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015671/R	MSN-D00CC000510	MISANO - Barra Cr-Co - d 5.5 mm L.510
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015672/R	MSN-D00CC000520	MISANO - Barra Cr-Co - d 5.5 mm L.520
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015673/R	MSN-D00CC000530	MISANO - Barra Cr-Co - d 5.5 mm L.530
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015674/R	MSN-D00CC000540	MISANO - Barra Cr-Co - d 5.5 mm L.540
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015675/R	MSN-D00CC000055	MISANO - Barra Cr-Co - d 5.5 mm L.55
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015676/R	MSN-D00CC000550	MISANO - Barra Cr-Co - d 5.5 mm L.550
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015677/R	MSN-D00CC000060	MISANO - Barra Cr-Co - d 5.5 mm L.60
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015678/R	MSN-D00CC000065	MISANO - Barra Cr-Co - d 5.5 mm L.65
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015679/R	MSN-D00CC000070	MISANO - Barra Cr-Co - d 5.5 mm L.70
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015680/R	MSN-D00CC000075	MISANO - Barra Cr-Co - d 5.5 mm L.75
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015681/R	MSN-D00CC000080	MISANO - Barra Cr-Co - d 5.5 mm L.80
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015683/R	MSN-D00CC000085	MISANO - Barra Cr-Co - d 5.5 mm L.85
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015684/R	MSN-D00CC000090	MISANO - Barra Cr-Co - d 5.5 mm L.90
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015685/R	MSN-D00CC000095	MISANO - Barra Cr-Co - d 5.5 mm L.95
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015686/R	MSN-D10CC00100	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.100
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015687/R	MSN-D10CC00110	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.110

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015688/R	MSN-D10CC00120	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.120
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015689/R	MSN-D10CC00130	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.130
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015690/R	MSN-D10CC00140	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.140
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015691/R	MSN-D10CC00150	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.150
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015692/R	MSN-D10CC00160	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.160
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015693/R	MSN-D10CC00170	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.170
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015694/R	MSN-D10CC00180	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.180
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015695/R	MSN-D10CC00190	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.190
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015696/R	MSN-D10CC00200	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.200
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015697/R	MSN-D10CC00210	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.210
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015698/R	MSN-D10CC00220	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.220
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015699/R	MSN-D10CC00230	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.230
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015700/R	MSN-D10CC00240	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.240
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015701/R	MSN-D10CC00025	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.25
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015702/R	MSN-D10CC00250	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.250
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015703/R	MSN-D10CC00030	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.30
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015704/R	MSN-D10CC00035	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.35
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015705/R	MSN-D10CC00040	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.40

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015706/R	MSN-D10CC00045	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.45
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015707/R	MSN-D10CC00050	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.50
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015708/R	MSN-D10CC00055	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.55
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015709/R	MSN-D10CC00060	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.60
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015710/R	MSN-D10CC00065	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.65
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015711/R	MSN-D10CC00070	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.70
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015712/R	MSN-D10CC00075	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.75
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015713/R	MSN-D10CC00080	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.80
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015714/R	MSN-D10CC00085	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.85
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015715/R	MSN-D10CC00090	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.90
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015716/R	MSN-D10CC00095	MISANO - Barra Cr-Co Precurvata - d 5.5 mm L.95
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015717/R	MSN-D10T500100	MISANO - Barra Precurvata - d 5.5 mm L.100
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015718/R	MSN-D10T500110	MISANO - Barra Precurvata - d 5.5 mm L.110
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015719/R	MSN-D10T500120	MISANO - Barra Precurvata - d 5.5 mm L.120
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015720/R	MSN-D10T500130	MISANO - Barra Precurvata - d 5.5 mm L.130
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015721/R	MSN-D10T500140	MISANO - Barra Precurvata - d 5.5 mm L.140
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015722/R	MSN-D10T500150	MISANO - Barra Precurvata - d 5.5 mm L.150
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015723/R	MSN-D10T500160	MISANO - Barra Precurvata - d 5.5 mm L.160

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015724/R	MSN-D10T500170	MISANO - Barra Precurvata - d 5.5 mm L.170
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015725/R	MSN-D10T500180	MISANO - Barra Precurvata - d 5.5 mm L.180
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015726/R	MSN-D10T500190	MISANO - Barra Precurvata - d 5.5 mm L.190
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015727/R	MSN-D10T500200	MISANO - Barra Precurvata - d 5.5 mm L.200
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015728/R	MSN-D10T500210	MISANO - Barra Precurvata - d 5.5 mm L.210
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015729/R	MSN-D10T500220	MISANO - Barra Precurvata - d 5.5 mm L.220
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015730/R	MSN-D10T500230	MISANO - Barra Precurvata - d 5.5 mm L.230
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015731/R	MSN-D10T500240	MISANO - Barra Precurvata - d 5.5 mm L.240
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015732/R	MSN-D10T500025	MISANO - Barra Precurvata - d 5.5 mm L.25
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015733/R	MSN-D10T500250	MISANO - Barra Precurvata - d 5.5 mm L.250
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015734/R	MSN-D10T500030	MISANO - Barra Precurvata - d 5.5 mm L.30
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015735/R	MSN-D10T500035	MISANO - Barra Precurvata - d 5.5 mm L.35
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015736/R	MSN-D10T500040	MISANO - Barra Precurvata - d 5.5 mm L.40
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015737/R	MSN-D10T500045	MISANO - Barra Precurvata - d 5.5 mm L.45
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015738/R	MSN-D10T500050	MISANO - Barra Precurvata - d 5.5 mm L.50
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015739/R	MSN-D10T500055	MISANO - Barra Precurvata - d 5.5 mm L.55
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015740/R	MSN-D10T500060	MISANO - Barra Precurvata - d 5.5 mm L.60
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015741/R	MSN-D10T500065	MISANO - Barra Precurvata - d 5.5 mm L.65

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015742/R	MSN-D10T500070	MISANO - Barra Precurvata - d 5.5 mm L.70
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015743/R	MSN-D10T500075	MISANO - Barra Precurvata - d 5.5 mm L.75
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015744/R	MSN-D10T500080	MISANO - Barra Precurvata - d 5.5 mm L.80
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015745/R	MSN-D10T500085	MISANO - Barra Precurvata - d 5.5 mm L.85
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015746/R	MSN-D10T500090	MISANO - Barra Precurvata - d 5.5 mm L.90
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015747/R	MSN-D10T500095	MISANO - Barra Precurvata - d 5.5 mm L.95
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015748/R	MSN-D00T500100	MISANO - Barra - d 5.5 mm L.100
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015749/R	MSN-D00T500105	MISANO - Barra - d 5.5 mm L.105
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015750/R	MSN-D00T500110	MISANO - Barra - d 5.5 mm L.110
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015751/R	MSN-D00T500115	MISANO - Barra - d 5.5 mm L.115
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015752/R	MSN-D00T500120	MISANO - Barra - d 5.5 mm L.120
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015753/R	MSN-D00T500125	MISANO - Barra - d 5.5 mm L.125
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015754/R	MSN-D00T500130	MISANO - Barra - d 5.5 mm L.130
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015755/R	MSN-D00T500135	MISANO - Barra - d 5.5 mm L.135
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015756/R	MSN-D00T500140	MISANO - Barra - d 5.5 mm L.140
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015757/R	MSN-D00T500145	MISANO - Barra - d 5.5 mm L.145
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015758/R	MSN-D00T500150	MISANO - Barra - d 5.5 mm L.150
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015759/R	MSN-D00T500155	MISANO - Barra - d 5.5 mm L.155

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015760/R	MSN-D00T500160	MISANO - Barra - d 5.5 mm L.160
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015761/R	MSN-D00T500165	MISANO - Barra - d 5.5 mm L.165
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015762/R	MSN-D00T500170	MISANO - Barra - d 5.5 mm L.170
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015763/R	MSN-D00T500175	MISANO - Barra - d 5.5 mm L.175
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015764/R	MSN-D00T500180	MISANO - Barra - d 5.5 mm L.180
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015765/R	MSN-D00T500185	MISANO - Barra - d 5.5 mm L.185
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015766/R	MSN-D00T500190	MISANO - Barra - d 5.5 mm L.190
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015767/R	MSN-D00T500195	MISANO - Barra - d 5.5 mm L.195
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015768/R	MSN-D00T500200	MISANO - Barra - d 5.5 mm L.200
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015769/R	MSN-D00T500210	MISANO - Barra - d 5.5 mm L.210
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015770/R	MSN-D00T500220	MISANO - Barra - d 5.5 mm L.220
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015771/R	MSN-D00T500230	MISANO - Barra - d 5.5 mm L.230
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015772/R	MSN-D00T500240	MISANO - Barra - d 5.5 mm L.240
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015773/R	MSN-D00T500025	MISANO - Barra - d 5.5 mm L.25
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015774/R	MSN-D00T500250	MISANO - Barra - d 5.5 mm L.250
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015775/R	MSN-D00T500260	MISANO - Barra - d 5.5 mm L.260
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015776/R	MSN-D00T500270	MISANO - Barra - d 5.5 mm L.270
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015777/R	MSN-D00T500280	MISANO - Barra - d 5.5 mm L.280

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015778/R	MSN-D00T500290	MISANO - Barra - d 5.5 mm L.290
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015779/R	MSN-D00T500030	MISANO - Barra - d 5.5 mm L.30
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015780/R	MSN-D00T500300	MISANO - Barra - d 5.5 mm L.300
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015781/R	MSN-D00T500310	MISANO - Barra - d 5.5 mm L.310
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015782/R	MSN-D00T500320	MISANO - Barra - d 5.5 mm L.320
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015783/R	MSN-D00T500330	MISANO - Barra - d 5.5 mm L.330
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015784/R	MSN-D00T500340	MISANO - Barra - d 5.5 mm L.340
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015785/R	MSN-D00T500035	MISANO - Barra - d 5.5 mm L.35
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015786/R	MSN-D00T500350	MISANO - Barra - d 5.5 mm L.350
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015787/R	MSN-D00T500360	MISANO - Barra - d 5.5 mm L.360
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015788/R	MSN-D00T500370	MISANO - Barra - d 5.5 mm L.370
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015789/R	MSN-D00T500380	MISANO - Barra - d 5.5 mm L.380
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015790/R	MSN-D00T500390	MISANO - Barra - d 5.5 mm L.390
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015791/R	MSN-D00T500040	MISANO - Barra - d 5.5 mm L.40
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015792/R	MSN-D00T500400	MISANO - Barra - d 5.5 mm L.400
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015793/R	MSN-D00T500410	MISANO - Barra - d 5.5 mm L.410
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015794/R	MSN-D00T500420	MISANO - Barra - d 5.5 mm L.420
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015795/R	MSN-D00T500430	MISANO - Barra - d 5.5 mm L.430

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015796/R	MSN-D00T500440	MISANO - Barra - d 5.5 mm L.290
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015797/R	MSN-D00T500045	MISANO - Barra - d 5.5 mm L.30
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015798/R	MSN-D00T500450	MISANO - Barra - d 5.5 mm L.300
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015799/R	MSN-D00T500460	MISANO - Barra - d 5.5 mm L.310
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015800/R	MSN-D00T500470	MISANO - Barra - d 5.5 mm L.320
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015801/R	MSN-D00T500480	MISANO - Barra - d 5.5 mm L.330
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015802/R	MSN-D00T500490	MISANO - Barra - d 5.5 mm L.340
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015803/R	MSN-D00T500050	MISANO - Barra - d 5.5 mm L.35
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015804/R	MSN-D00T500500	MISANO - Barra - d 5.5 mm L.350
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015805/R	MSN-D00T500510	MISANO - Barra - d 5.5 mm L.360
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015806/R	MSN-D00T500520	MISANO - Barra - d 5.5 mm L.370
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015807/R	MSN-D00T500530	MISANO - Barra - d 5.5 mm L.380
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015808/R	MSN-D00T500540	MISANO - Barra - d 5.5 mm L.390
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015809/R	MSN-D00T500055	MISANO - Barra - d 5.5 mm L.40
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015810/R	MSN-D00T500550	MISANO - Barra - d 5.5 mm L.400
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015811/R	MSN-D00T500060	MISANO - Barra - d 5.5 mm L.410
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015812/R	MSN-D00T500065	MISANO - Barra - d 5.5 mm L.420
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015813/R	MSN-D00T500070	MISANO - Barra - d 5.5 mm L.430

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015814/R	MSN-D00T500075	MISANO - Barra - d 5.5 mm L.75
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015815/R	MSN-D00T500080	MISANO - Barra - d 5.5 mm L.80
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015816/R	MSN-D00T500085	MISANO - Barra - d 5.5 mm L.85
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015817/R	MSN-D00T500090	MISANO - Barra - d 5.5 mm L.90
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015818/R	MSN-D00T500095	MISANO - Barra - d 5.5 mm L.95
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015819/R	MSN-E00T500011	MISANO - Uncino Corpo Allungato - Taglia 11 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015820/R	MSN-E00T500007	MISANO - Uncino Corpo Allungato - Taglia 7 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015821/R	MSN-E00T500009	MISANO - Uncino Corpo Allungato - Taglia 9 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015822/R	MSN-E10T500011	MISANO - Lumbar angled hook size 11 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015823/R	MSN-E10T500007	MISANO - Lumbar angled hook size 7 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015824/R	MSN-E10T500009	MISANO - Lumbar angled hook size 9 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015825/R	MSN-E20T500011	MISANO - Lumbar narrow hook size 11 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015826/R	MSN-E20T500007	MISANO - Lumbar narrow hook size 7 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015827/R	MSN-E20T500009	MISANO - Lumbar narrow hook size 9 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015828/R	MSN-E30T500011	MISANO - Lumbar wide hook size 11 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015829/R	MSN-E30T500007	MISANO - Lumbar wide hook size 7 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015830/R	MSN-E30T500009	MISANO - Lumbar wide hook size 9 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015831/R	MSN-E40T510007	MISANO - Oblique hook size 7 DX con Dado

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015832/R	MSN-E40T520007	MISANO - Oblique hook size 7 SX con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015833/R	MSN-E50T510008	MISANO - Offset hook DX con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015834/R	MSN-E50T520008	MISANO - Offset hook SX con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015835/R	MSN-E60T500005	MISANO - Pedicle hook size 5 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015836/R	MSN-E60T500007	MISANO - Pedicle hook size 7 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015837/R	MSN-E60T500009	MISANO - Pedicle hook size 9 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015838/R	MSN-E70T500005	MISANO - Thoracic laminar hook size 5 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015839/R	MSN-E70T500007	MISANO - Thoracic laminar hook size 7 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015840/R	MSN-E70T500009	MISANO - Thoracic laminar hook size 9 con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015841/R	MSN-F10T500020	MISANO - Angled offset L20 mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015842/R	MSN-F10T500030	MISANO - Angled offset L30 mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015843/R	MSN-F10T500060	MISANO - Angled offset L60 mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015844/R	MSN-F00T500020	MISANO - Offset L20 mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015845/R	MSN-F00T500030	MISANO - Offset L30 mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015846/R	MSN-F00T500060	MISANO - Offset L60 mm con Dado
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015847/R	MSN-G00T300000	MISANO - Dominò con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015848/R	MSN-G10T300000	MISANO - Elemento di connessione longitudinale con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015849/R	MSN-H00T305515	MISANO - Cross-link rigido piatto d 5,5 15 mm con Grani

CND	RDM	CODE	DESCRIPTION
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015850/R	MSN-H00T305518	MISANO - Cross-link rigido piatto d 5,5 18 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015851/R	MSN-H00T305521	MISANO - Cross-link rigido piatto d 5,5 21 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015852/R	MSN-H00T305524	MISANO - Cross-link rigido piatto d 5,5 24 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015853/R	MSN-H00T305527	MISANO - Cross-link rigido piatto d 5,5 27 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015854/R	MSN-H00T305530	MISANO - Cross-link rigido piatto d 5,5 30 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015855/R	MSN-H00T305534	MISANO - Cross-link rigido piatto d 5,5 33 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015856/R	MSN-H00T305536	MISANO - Cross-link rigido piatto d 5,5 36 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015857/R	MSN-H00T305539	MISANO - Cross-link rigido piatto d 5,5 39 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015858/R	MSN-H00T305542	MISANO - Cross-link rigido piatto d 5,5 42 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015859/R	MSN-H00T305545	MISANO - Cross-link rigido piatto d 5,5 45 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015860/R	MSN-H00T305548	MISANO - Cross-link rigido piatto d 5,5 48 mm con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015861/R	MSN-I00T300000	MISANO - Cross-link modulare - short con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015862/R	MSN-I00T300001	MISANO - Cross-link modulare - medium con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015863/R	MSN-I00T300002	MISANO - Cross-link modulare - long con Grani
P09070302 - SYSTEMS FOR DORSO-LUMB-SACRAL FIXATION	2015864/R	MSN-I00T300003	MISANO - Cross-link modulare - extra long con Grani